



Intellihot

www.intellihot.com
2900 W. Main Street, Galesburg IL 61401
877-835-1705
support@intellihot.com

January 29, 2021
Bulletin#: TB2021-01

Technical Bulletin regarding strainer screen change - Neurons

Models: iN401 & iN501

Description:

Effective with all Neurons manufactured after January 26, 2021, the wire mesh in the supplied inlet water strainers for the Neurons (models iN401 & iN501) is being upgraded. The previous supplied strainer included a No. 20 mesh (which allowed particles up to 841 microns to pass through), and the upgraded supplied strainer will include a No. 30 mesh (allowing particles up to 595 microns).

While the mesh upgrade will catch more particulate, this new design still provides the exact same amount of free opening and thus will not affect the performance or flow in any way.

Part number for new strainer is: IGT-FTT0282.

To contact Intellihot product support, please call us at 877-835-1705 or email us at support@intellihot.com